

## **SERVOENT neo NICU Ventilator**

1. Easiest possible setting
  2. Can be operated by without much knowledge of ventilator
  3. Fastest Ready to ventilate
  4. NIV mode inbuilt
  5. Bubble CPAP inbuilt
  6. It can be portable for Ambulance use
  7. Can be run on any one gas in emergency
  8. Battery backup
  9. One button ventilation .....
- ....many more for your baby



### **MEDI ELECTRONICS**

G-20, J-2 Building, Shirdi Nagar, Bhayander (E),  
Dist. Thane, Mumbai - 401 105, Maharashtra, India.  
B: 022 28186268 • E: [info@medielectronics.in](mailto:info@medielectronics.in)  
W: [www.medielectronics.in](http://www.medielectronics.in)



## **MEDI ELECTRONICS**

*...."Doctor's convenience, patient's comfort!"*



## **EQUIPMENTS YOU CAN RELY ON**

## DIGITAL TOURNIQUET

Model: T1



- ♦ Light Weight 1.3Kg
- ♦ Battery Backup 5 Hrs
- ♦ Compact
- ♦ User friendly white LCD screen
- ♦ Silicon Bladder cuff
- ♦ Fully Microprocessor controlled
- ♦ Service & Calibration from front panel key
- ♦ Timer can also individually used
- ♦ Set and Actual both display continually
- ♦ Easy to use
- ♦ Settings can be changed while on going operation
- ♦ Settings can be done in stand by position also
- ♦ Works on mains as well battery also
- ♦ No valve noise
- ♦ User can enable or disable. Check cuff alarm
- ♦ Color coded cuff set
- ♦ Silicon extension tube



## AIR COMPRESSOR FOR VENTILATOR

Model: P60SIDigital

- ♦ Microprocessor Controlled
- ♦ Diagnostic display and alarming for MCB tripper, All fans, Temperature, Water draining etc.
- ♦ Auto Drain Inbuilt
- ♦ Only 30 kg Max Weight
- ♦ Compact size
- ♦ Thomas Motor used
- ♦ Noise 45db at 1mtr distance
- ♦ 4 Bar ( 60psi ) regulated output



## OXYGEN MONITOR

Model: O2 100

- ♦ Fast response
- ♦ Accurate & Reliable
- ♦ Auto Calibration
- ♦ Mains Operated
- ♦ 8 hrs. Battery Backup
- ♦ User selectable audiovisual alarm limits
- ♦ Large LCD Display
- ♦ Size 17 x 6 x 13



## CUFF SET FOR Digital Tourniquet

- ♦ Superior fabric
- ♦ Silicon Badder
- ♦ Color Coded
- ♦ Long Lasting
- ♦ Reusable



## OXYGEN PIPELINE PRESSURE MONITOR

Model: Gasm100



- ♦ Fast Response
- ♦ Accurate & Reliable
- ♦ Auto calibration
- ♦ Operated on main supply
- ♦ User selectable audiovisual alarm limits
- ♦ Both Upper & Lower alarm limits selectable by only two key
- ♦ Large Backlight LCD
- ♦ Compact size
- ♦ Wall Mountable
- ♦ Battery Backup(optional)



## PROXIMAL AIRWAY PRESSURE MONITOR

Model: PIP100

- ♦ Digit as well Bar Graphic Display
- ♦ User selectable Audiovisual alarm limits
- ♦ With Active PIP Relief Valve
- ♦ Accurate & Reliable
- ♦ Auto Calibration
- ♦ Very Fast Reading
- ♦ Operated on main supply
- ♦ 8 hrs Battery Backup(optional)
- ♦ Both Upper & Lower alarm limits selectable by only two key
- ♦ Size in cm 17 x 6 x 13



## OXYGEN PIPELINE PRESSURE MONITOR

Model: ME-150

- ♦ Fast Response
- ♦ Accurate & Reliable
- ♦ Max pressure 150 psi
- ♦ Operated on main supply
- ♦ ABS Plastic body
- ♦ Easy to read pressure gauge dial
- ♦ Compact size
- ♦ Wall Mountable



## DIGITAL TIMER SWITCH

Model: ME10A

- ♦ LCD Interface
- ♦ Accurate & Reliable
- ♦ Easy Setup
- ♦ Compact size
- ♦ Wall Mountable
- ♦ Up to 10amp switching
- ♦ ABS Plastic body

