


Proximal Airway Pressure Monitor



FEATURES:

- Digit as well Bar graphic display.
- User selectable audiovisual alarm limits.
- With active PIP relief valve.
- Accurate & reliable.
- Auto calibration.
- Very fast reading.
- Operated on main supply.
- 8 hrs. Battery backup (optional).
- Both upper & lower alarm limits selectable by only two key.
- Size in cm 17 x 6 x 13.
- Intended to use.

Technical Specification:

Ambient operating temperature Range:	10°C to 50°C.
Operating Altitude:	-557 ft. to +16,000 ft / -170m to +4877m.
Patient pressure Range:	0 to 100 cmH ₂ O.
Sample Pressure inlet port size:	5 mm male connector.
Transport/Storage Temperature:	-13°F to +131°F (-25°C to +55°C).
Accuracy:	+/- 1 cmH ₂ O.
Humidity:	15% to 95% non-condensing.

Dimensions:	(17cm x 6cm x 13cm).
Weight:	0.85 Kg.
Power:	230V AC (10W)
Equipment Classification:	Class 2 Device, Continuous Operation
Water protection:	Not protected
Alarm:	Audio Visual
Alarm range:	default low = 0, High = 35. User selectable low = 0 TO 20. User selectable High = 5 TO 90.
Display:	LCD.
Battery Backup (Optional):	8 Hrs.